**ENGLISH** SPECIFICATION

1 test/kit, 5 tests/kit, 7 tests/kit, 10 tests/kit, 15 tests/kit,  
20 tests/kit, 25 tests/kit

CE 1434 For self-testing

IVD In vitro diagnostic medical device

**INTENDED USE**

The GLINE-2019-nCoV Ag Test is a colloidal gold immunoassay (CGIA) for qualitative detection of SARS-CoV-2 nucleocapsid antigens in direct nasal swab from individuals who are suspected of SARS-CoV-2 within the first seven days of the onset of symptoms. The GLINE-2019-nCoV Ag Test does not differentiate between SARS-CoV and SARS-CoV-2.

Positive results indicate that viral antigens are present in the sample. Positive results do not exclude the possibility of bacterial infection or co-infection with other viruses.

Negative results do not exclude SARS-CoV-2 infection. SARS-CoV-2 infection or not should at least take into account your disease history clinical signs and symptoms consistent with COVID-19.

The GLINE-2019-nCoV Ag is intended for self-testing. Teenagers and children under 18 years old should follow the guidance from their guardians to perform a self-test. Elders above 80 should seek for assistance to perform a self-test.

**SUMMARY AND EXPLANATION**

The novel coronaviruses belong to the β genus. COVID-19 is an acute respiratory infectious disease. People are generally susceptible. Currently, the patients infected by the novel coronavirus are the main source of infection; asymptomatic infected people can also be an infectious source. Based on the current epidemiological investigation, the incubation period is 1 to 14 days, mostly 3 to 7 days but still remains infectious during this period. The main manifestations include fever, fatigue and dry cough. Nasal congestion, runny nose, sore throat, myalgia and diarrhea are found in a few cases.

**PRINCIPLE OF THE TEST**

The GLINE-2019-nCoV Ag Test employs colloidal gold immunoassay technology in a sandwich design to detect the SARS-CoV-2 nucleocapsid protein.

- When specimens are processed and added into sample well, the specimen will flow along the test strip by capillary action. If the SARS-CoV-2 nucleocapsid protein exists in the specimen, the protein will combine with colloidal gold labeled SARS-CoV-2 antibodies. This immune complex will be captured by immobilized anti-SARS-CoV-2 antibodies in the detection region (T) and form a coloured detection line, which means SARS-CoV-2 antigen positive.
- The detection cassette also includes a control region (C). Upon completion of a valid test, a coloured control line should appear regardless of whether the SARS-CoV-2 antigen is present in the sample. If the control line does not appear, the test result is invalid and the sample should be retested with another test cassette.

**KIT COMPONENTS**

GLINE-2019-nCoV Ag contains:

REF	G86 256	G86 257	G86 278	G86 279	G86 280	G86 273	G86 281
Test Cassette (pcs)	1	5	7	10	15	20	25
Extraction Buffer, 400 µL/tube (pcs)	1	5	7	10	15	20	25
Sterile Swabs (pcs)	1	5	7	10	15	20	25
Biohazard Sample Bags (pcs)	1	5	7	10	15	20	25
Instruction For Use (pcs)	1	1	1	1	1	1	1

The test cassette is packed in a sealed foil pouch with 1 Desiccant Pouch and composed with (see FIG. 1):

- Sample pad: glass fiber membrane.
- Gold pad: colloidal gold particles labeled mouse anti-SARS-CoV-2 antibodies and peptides coupled with colloidal gold.
- Nitrocellulose membrane: nitrocellulose (NC) membrane with immobilized mouse anti-SARS-CoV-2 antibodies (line T) and mouse anti-peptide antibodies (line C).
- Absorbent pad: cellulose-acetate membrane.
- Backing pad: PVC card.
- Housing: PS (polystyrene) plastic.


FIG. 1

**MATERIALS REQUIRED (BUT NOT PROVIDED)**

- Timer.

**WARNINGS AND PRECAUTIONS**

IVD For in vitro diagnostic use only.

- This test has been authorized only for the detection of proteins from SARS-CoV-2, not for any other virus or pathogen.
- Use immediately after opening the pouch containing the test cassette.
- Avoid touching any bleeding areas of the nostril area during specimens collection, as excess blood or mucus on the swab may interfere with the test and yield a false result.
- Do not use if the test device package is damaged.
- Do not use the kit contents beyond the expiration date.
- Do not eat, drink, or smoke in the area where the specimens and kit contents are handled.
- Do not interchange kit contents from different lots.
- Do not make any decision of medical relevance without first consulting your medical practitioner.
- Clean and disinfect all spills of specimens or buffer using a suitable disinfectant.
- If the extraction buffer makes contact with the skin or eye, wash/ flush with a large volume of water. If skin irritation, rash or other abnormal reaction occurs, please get medical advice/attention.
- Follow the instruction strictly when performing a self-test. Children and elders please use the test under the guardian assistance.

**REAGENT HANDLING**

- The IFU should be followed precisely for accurate result.
- Discard after first use. The test cannot be used more than once.
- Do not touch the test area of test strip.

**STORAGE AND STABILITY**


- The product should be stored at 2-30°C, dry and out of light, its validity for 18 months.
- Each test cassette should be used within 1 hour of opening the sealed foil pouch. The test should be used in a room temperature (20-30°C) environment with a humidity of no more than 75%.
- If the test has been stored in a refrigerator, allow the test to equilibrate to room temperature for 30 minutes before use.
- The test should be used immediately after being removed from its sealed foil pouch.

**SPECIMEN COLLECTION AND HANDLING**

- It is essential that correct specimen collection and preparation methods should be followed. Improper collection or handling of samples may yield false results.

• **Specimen Collection:**

- To collect a nasal swab sample, carefully insert the swab into the nostril exhibiting the most visible drainage, or the nostril that is most congested if drainage is not visible. Using gentle rotation, push the swab (about 1.5cm into the nostril). Rotate the swab 5 times or more against the nasal wall then slowly remove from the nostril. Using the same swab, repeat sample collection in the other nostril.

• **Specimen Transport and Storage:**


- After swabbing, process the swab in the extraction buffer as soon as possible. Do not place the swab back into the swab packaging sleeve after specimen collection.
- Specimens should be tested within 30 minutes. Do not freeze or transport the sample for later testing.
- To ensure optimal test performance and avoid possible contamination, collect the specimen and finish the test one after one or perform the test under social distance if there are multiple members in the same facility. All specimens must be handled with caution to prevent cross contamination.

**ASSAY PROCEDURE**

Read the IFU carefully before running a test.

Step 1:

- Peel off the foil film of the extraction tube which contains buffer. Set the extraction tube on the stand hole of the package box.

Step 2:

- After sample collection, insert the swab into the sample extraction tube, plunge and press the swab tip up and down in the fluid for a minimum of 15 seconds, taking care not to splash contents out of the tube.

Step 3:


- Remove the swab while squeezing the sides of the tube to extract the liquid from the swab.

Step 4:

- Press the included dispensing tip onto the extraction tube containing the processed sample. Mix thoroughly by swirling or flicking the bottom of the tube.

Step 5: Adding the specimen to the test cassette

- Invert the extraction tube and hold it vertically.
- Gently squeeze the ridged body of the tube, dispensing three (3) drops of the processed specimen into the sample well.

Step 6:

- Read the results which should be ready at 15-20 minutes after adding the processed sample to the sample well. Results read after 20 minutes are invalid.


Step 7:

- After reading the results, place all the components of the test kit in the provided biohazard sample bag and dispose according to local regulation. Do not reuse any used components.

## INTERPRETATION OF RESULTS


### Positive (+)

- A coloured line appear in the control region (C) and the detection region (T). This indicates that the SARS-CoV-2 antigens are present in the sample. The detection region (T) can be very faint. Any coloured line visible here is positive.
- If you get a positive result with the GLINE-2019-nCoV Ag Test, it's very likely you have COVID-19. You should follow local guidelines for self-isolation and immediately contact your general practitioner/doctor or the local health department in accordance with the instructions of your local authorities.


### Negative (-)

- A coloured line appears in the control region (C) but not the detection region (T). This indicates that the SARS-CoV-2 antigens are not present in the sample.
- You are unlikely to have COVID-19. However, a negative result does not rule out COVID-19. This means you could possibly still have COVID-19 even though the test result is negative. If you experience symptoms such as headaches, migraines, fever, loss of sense of smell or taste, contact the nearest medical facility according to the rules of your local authority. In addition, you can repeat the test with a new test kit. In case of suspicion, repeat the test after 1-2 days, as the coronavirus cannot be precisely detected in all phases of an infection. In consideration of the pandemic, it is best to self-isolate before the next test.

### Invalid

- No coloured line appears in the control region (C) after performing the test. The directions may not have been followed correctly or the test may have deteriorated. You need review the instruction for use again and repeat the test with a new cassette.

## LIMITATIONS

- This test has been evaluated for use with human specimen material only.
- Positive test results from this test do not differentiate between SARS-CoV and SARS-CoV-2.
- The accuracy of the test relies on the whole testing procedure being performed correctly. Failure to follow the instructions will adversely affect test accuracy. Retesting is recommended to confirm the results if the test results are inconsistent with clinical symptoms.
- This product is only used for qualitative detection of SARS-CoV-2 antigens in nasal swab specimens and cannot detect antigen concentration in samples.
- Results from this are for clinical reference only, and should not be used as the sole basis for clinical diagnosis and treatment.

- Please take the necessary safety measures (e.g. face mask, gloves) when testing for other people.
- Wash hands thoroughly or use hand sanitizer after handling.
- Wash or sanitize your hands and dry them before sampling and testing
- Negative results in cases where COVID-19 disease is suspected should be reviewed by a confirmatory test if necessary.
- Possibility of false negative result may be:
  - 1) Improper collection, handling, transport of samples, or low virus in samples may also cause false negative results.
  - 2) The sample may have been collected when the patient's virus is low, or the virus in the region from which the sample was collected is not enough.

## PERFORMANCE CHARACTERISTICS

### Clinical performance

A limited study involving three hundred and sixty-five (365) direct nasal swabs was performed. The samples were sequentially enrolled from two locations and tested immediately after collection. Results from the GLINE-2019-nCoV Ag Test were compared with a reference RT-PCR assay.

GLINE-2019-nCoV Ag Results	Reference PCR Results		
	Positive	Negative	Total
Positive	133	1	134
Negative	5	226	231
Total	138	227	365

Sensitivity of GLINE-2019-nCoV Ag Test: 133/138= 96.38%(95%CI:91.75% - 98.81%) . In total 138 PCR confirmed positive samples: 133 PCR confirmed positive samples were correctly detected by COVID-19 Antigen Rapid Test. There are 5 false negative cases.

Specificity of GLINE-2019-nCoV Ag Test: 226/227= 99.56% (95%CI:97.57% - 99.99%). In total 227 PCR confirmed negative samples: 226 PCR confirmed negative samples were correctly detected by COVID-19 Antigen Rapid Test. There are 1 false positive cases.

### Limit of Detection

The limit of Detection (LoD) of the GLINE-2019-nCoV Ag Test was determined to be the lowest concentration resulting in positive detection of twenty (20) out of twenty (20) replicates. Based on this testing the LoD was confirmed as:  $2.0 \times 10^2$  TCID<sub>50</sub>/mL.

### Cross-Reactivity

Cross-reactivity of the GLINE-2019-nCoV Ag Test was evaluated by testing a panel of high prevalence respiratory pathogens that could potentially cross-react with the GLINE-2019-nCoV Ag Test. The following pathogens were not detected using GLINE-2019-nCoV Ag Test.

Potential Cross-Reactant	Concentration Tested
Human coronavirus HKU1	$1.0 \times 10^5$ TCID <sub>50</sub> /mL
Human coronavirus OC43	$1.0 \times 10^5$ TCID <sub>50</sub> /mL
Human coronavirus 229E	$1.0 \times 10^5$ TCID <sub>50</sub> /mL
Human coronavirus NL63	$1.0 \times 10^5$ TCID <sub>50</sub> /mL
Influenza A H1N1	$2.0 \times 10^5$ TCID <sub>50</sub> /mL
Influenza A H3N2	$2.5 \times 10^5$ TCID <sub>50</sub> /mL
Influenza B	$3.0 \times 10^5$ TCID <sub>50</sub> /mL

Respiratory syncytial virus	4.0 x 10 <sup>5</sup> TCID <sub>50</sub> /mL
Rhinovirus	2.0 x 10 <sup>5</sup> PFU/mL
Adenovirus	1.0 x 10 <sup>5</sup> TCID <sub>50</sub> /mL
Enterovirus	4.0 x 10 <sup>5</sup> TCID <sub>50</sub> /mL
Parainfluenza	1.0 x 10 <sup>5</sup> TCID <sub>50</sub> /mL
Mycoplasma pneumoniae	1.0 x 10 <sup>6</sup> CFU/mL
MERS coronavirus	1.0 x 10 <sup>5</sup> TCID <sub>50</sub> /mL

## Anti-interference

A study was performed demonstrate that eighteen potentially interfering substances that may be found in the upper respiratory tract do not cross-react or interfere with the detection of SARS-CoV-2 in the GLINE-2019-nCoV Ag Test.

Interfering Substance	Concentration Tested
Mucin	5 mg/mL
Blood (human)	5% v/v
Zanamivir	0.1 mg/mL
Oseltamivir	2.5 µg/mL
Tobramycin	1.0 mg/mL
Interferon alpha	3 x 10 <sup>5</sup> IU/mL
Ribavirin	0.4 mg/mL
Paramivir	0.15 mg/mL
Lopinavir	0.2 mg/mL
Levofloxacin	0.2 mg/mL
Azithromycin	1 µg/mL
Ritonavir	0.3 mg/mL
Histamine hydrochloride	0.1mg/mL
Sodium Chloride	5% v/v
Oxymetazoline	15% v/v
Dexamethasone	1.2 mg/dL
Triamcinolone	5% v/v
Fluticasone Propionate	5% v/v

## Hook Effect

No high dose hook effect was observed up to 1.6×10<sup>5</sup> TCID<sub>50</sub>/mL of heat inactivated SARS-CoV-2 virus with the GLINE-2019-nCoV Ag Test.

## REFERENCES

- Zhou P, Yang XL, Wang XG, et al. A pneumonia outbreak associated with a new coronavirus of probable bat origin[J]. Nature, 2020.
- Lu R, Zhao X, Li J, et al. Genomic characterisation and epidemiology of 2019 novel coronavirus: implications for virus origins and receptor binding[J]. Lancet, 2020.

SHENZHEN YHLO BIOTECH CO., LTD.  
Building 1, YHLO Biopark, Baolong 2nd Road,  
Baolong Subdistrict, Longgang District, 518116  
Shenzhen, PEOPLE'S REPUBLIC OF CHINA

Lotus NL B.V.  
Koningin Julianaplein 10, 1e Verd, 2595AA, The  
Hague, Netherlands.  
E-mail: peter@lotusnl.com

Distributed by BGI  
www.bgi.com/global/  
info@bgi.com

## ANNEX A:

### Explanation of abbreviation

Abbreviation	Explanation
	Catalogue number
	Contains sufficient for <n> tests
	Manufacturer
	Authorized representative in the European Community
	Caution
	Consult instructions for use
	<i>In vitro</i> diagnostic medical device
	Batch code
	Date of manufacture
	Use-by date
	Temperature limit (2-30°C)
	This way up
	Do not re-use
	Do not use if package is damaged
	Keep dry
	Keep away from sunlight
	CE Conformity Marking
	Distributor

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1 test/kit, 5 tests/kit, 7 tests/kit, 10 tests/kit, 15 tests/kit, 20 tests/kit, 25 tests/kit

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### REFERENCIAS

- Zhou P, Yang XL, Wang XG, et al. *A pneumonia outbreak associated with a new coronavirus of probable bat origin*[J]. Nature, 2020.
- Lu R, Zhao X, Li J, et al. *Genomic characterisation and epidemiology of 2019 novel coronavirus: implications for virus origins and receptor binding*[J]. Lancet, 2020.

 SHENZHEN YHLO BIOTECH CO., LTD.  
Building 1, YHLO Biopark, Baolong 2nd Road,  
Baolong Subdistrict, Longgang District, 518116  
Shenzhen, PEOPLE'S REPUBLIC OF CHINA

 Lotus NL B.V.  
Koningin Julianaplein 10, 1e Verd, 2595AA,  
The Hague, Netherlands.  
E-mail: peter@lotusnl.com

 Distributed by BGI  
www.bgi.com/global/  
info@bgi.com

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